



**Leader in providing highest quality RE-Furbished Ultrasound**

## Inventory

(December 15, 2008)  
(Rev 1)



**TOSHIBA**



### Shipping Guidelines

**Zone A – 3-5 Days – In Showroom**

**Zone B – 10 Days – In Warehouse**

**Zone C – 30 Days – In Transit**

Portable Systems – Page 2- 4  
Stand-Alone Systems – Page 4 - 14  
Printer Accessories – Page 15

## PORTABLE SYSTEMS

### ACUSON / SIEMENS

- 1 **Cypress PLUS – Rev. 20 w/ Stress Echo & DICOM (MFG: 2007) -- NEW in the BOX**  
Cardiac/Vascular calc., Cypress Viewer, DICOM Networking, Stress Echo  
Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG Cable, MOD, IMT Option  
**3V2C – Adult Cardiac Sector, 7L3 – Vascular Linear 1 year warranty**  
**ZONE A**
- 2 **Cypress PLUS – Rev. 20 w/ Stress Echo & DICOM (MFG: 2007)-- DEMO**  
Cardiac/Vascular calc., Cypress Viewer, DICOM Networking, Stress Echo  
Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG Cable, MOD, IMT Option  
**3V2C – Adult Cardiac Sector, 7L3 – Vascular Linear**  
**ZONE C**
- 3 **Cypress – Rev. 12.31 KPI 24825**  
Cardiac Calculation, DICOM Networking, Stress Echo  
Color Doppler, CW/PW Doppler, M-Mode, ECG Cable, MOD  
**3V2C – Adult Cardiac Sector**  
**ZONE A PENDING SALE**  
  
**Option: 4C1 Convex - (with system sale only)**  
**Cypress Deluxe Cart- reconditioned (For rev. 14 and below)**

### GENERAL ELECTRIC (GE)

- 4 **LOGIQ 100 Alpha - B/W – OB System - PORTABLE**  
**C36 (3.5MHz) Convex**  
**ZONE B**
- 5 **LOGIQ 100 MP - B/W – OB System (MFG: 2001) - PORTABLE**  
**C36 (3.5MHz) Convex**  
**ZONE A**
- 6 **LOGIQ 100 PRO - B/W – OB System (MFG: 2003) - PORTABLE**  
**C36 (3.5MHz) Convex**  
**ZONE A**
- 7 **LOGIQBOOK – OB/Vascular/Small Parts – Rev. 2.03 (MFG: 2003)**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics,  
**3C-RS (2-5MHz) Convex**  
**(110/220V – PAL/NTSC) – ZONE A**  
**“Special price” 2 units at this price**
- 8 **LOGIQBOOK – OB/Vascular/Small Parts – Rev. 2.03 (MFG: 2003/2004)**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics,  
**3C-RS (2-5MHz) Convex**  
**(110/220V – PAL/NTSC) – ZONE A**
- 9 **LOGIQBOOK – OB/Vascular/Small Parts – Rev. 2.03 (MFG: 2004)**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics,  
**3C-RS (2-5MHz) Convex, E8C-RS (4-9MHz) Endo-cavity OR 8L-RS Linear**  
**“NEW PROBES” (110/220V – PAL/NTSC) – ZONE A**
- 10 **LOGIQBOOK XP – Cardiac/OB/Vascular/Small Parts – Rev. 2.1 (MFG: 2004)**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics, DICOM  
**3C-RS (2-5MHz) Convex, E8C-RS (4-9MHz) Endo-cavity OR 8L-RS (4-10MHz) Linear**  
**“NEW PROBES” (110/220V – PAL/NTSC) – ZONE A**
- 11 **LOGIQBOOK XP – Cardiac/OB/Vascular/Small Parts – Rev. 2.1 (MFG: 2005)**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics, DICOM  
**3C-RS (2-5MHz) Convex, E8C-RS (4-9MHz) Endo-cavity OR 8L-RS (4-10MHz) Linear**  
**“NEW PROBES” (110/220V – PAL/NTSC) – ZONE A**

- 12 **LOGIQBOOK XP – Cardiac/OB/Vascular/Small Parts – Rev. 2.1 (MFG: 2007)-DEMO UNIT**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics, Report Writer  
With choice of (2) NEW transducers ... 3C-RS, 8L-RS, E8C-RS, 3S-RS (cardiac Sector)  
“NEW PROBES” (110/220V – PAL/NTSC) – ZONE C (2) Units available
- 13 **LOGIQBOOK XP VET – Veterinary Package – Rev. 2.1 (MFG: 2004)**  
Color Doppler, PW Doppler, Angio, M-Mode, ATO, Tissue Harmonics, DICOM  
8C-RS (4-9MHz) Micro Convex  
(110/220V – PAL/NTSC) ZONE A

**Additional Transducer:**

- 3C-RS (2-5MHz) Convex - new
  - 4C-RS (1.7-4.0MHz) Convex - new (for LOGIQ-e/ Logiq-i/Voluson-i)
  - E8C-RS (4-10MHz) Endo-cavity used / new
  - 8C-RS (4-10MHz) Micro-convex for Neonatal - new
  - 8L-RS (4-10MHz) Linear - used / new (for XP only)
  - 3S-RS (1.5-3.5MHz) Cardiac Sector - used/new (for XP only)
  - i739L-RS (4-10MHz) Veterinary/Intraoperative Linear - new (for XP only)
  - 12L-RS (5-12MHz) Linear - new – (for LOGIQ-e/ Logiq-i/Voluson-i)
- 14 **LOGIQ-e BTO6 (Portable) – CVO**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics, **Dicom**  
Speckle Reduction Imaging (SRI), Auto Tissue Optimization, OB/Vascular/Radiology Pkg.  
4C-RS (2-5MHz) Convex (new), E8C-RS (4-10MHz) Endo-cavity (new)  
**ZONE A**  
“Special price”
- 15 **LOGIQ-e BTO6 MFG: 2006-2007 (Portable) – CVO KPI 24888**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics, Easy 3D, **Dicom**  
Speckle Reduction Imaging (SRI), Auto Tissue Optimization, OB/Vascular/Radiology Pkg.  
4C-RS (2-5MHz) Convex (new)  
**ZONE A**
- 16 **LOGIQ-e BTO6 MFG: 2006-2007 (Portable) – CVO KPI 24889**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics, Easy 3D, **Dicom**, ECG  
Speckle Reduction Imaging (SRI), Auto Tissue Optimization, OB/Vascular/Radiology Pkg.  
4C-RS (2-5MHz) Convex (new), 12L-RS (5-12MHz) Linear (new), 3S-RS (1.5-3.5Mhz) Cardiac(new)  
**ZONE A**
- 17 **LOGIQ-e BTO6 MFG: 2007 (Portable) – CVO**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics, ECG, Speckle Reduction Imaging (SRI),  
Auto Tissue Optimization, OB/Vascular/Radiology Pkg. LOGIQVIEW, **DICOM**, Spatial Compounding,  
Anatomical M-Mode (AMM)  
4C-RS (2-5MHz) Convex (new)  
**ZONE C**
- 18 **VOLUSON-e BT06 MFG: 2006 (Portable) – Realtime 4D KPI 24914 / 24913**  
Color Doppler, PW Doppler, M-Mode, Angio, Tissue Harmonics, ATO  
DICOM 3.0, Realtime 4D, Multi-plane 3D, Surface 3D, OB/Vascular/Radiology Pkg.  
RAB2-5RS (2-5MHz) 4D Convex, RIC5-9RS (5-9MHz) 4D Endo-cavity  
**ZONE A**
- 19 **VIVID-i2 BT06 MFG: 2006-2007 (Portable) – Cardiac/Vascular Pkg. System**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics, DICOM Media, DVD-RW, Echo Pac, ECG  
Stress Echo, MPEG VUE, Extended Delivery  
3S-RS (1.5-3.5MHz) Cardiac Sector , 8L-RS (4-10MHz) Vascular Linear  
**ZONE A**
- 20 **VIVID-i BT06 MFG: 2006 (Portable) – CardioVascularPkg. System**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics,  
DICOM Media, DVD-RW, Echo Pac, ECG, **Stress Echo**, MPEG VUE  
3S-RS (1.5-3.5MHz) Cardiac Sector , 8L-RS (4-10MHz) Vascular Linear  
**ZONE A**

- 21 **VIVID-i BT06 MFG: 2006 (Portable) – CardioVascularPkg. System**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Tissue Harmonics,  
DICOM Media, DVD-RW, Echo Pac, ECG, **Stress Echo**, MPEG VUE  
**3S-RS (1.5-3.5MHz) Cardiac Sector , 7S (4-9MHz) Pediatric Cardiac Sector**  
**ZONE A**

### **MEDISON**

- 22 **SA-PICO Portable Digital Color Doppler (MFG: 2007)**  
Color Doppler, PW Doppler, M-Mode, OB/GYN/Vascular Pkg., SonoView, Surface 3D, DICOM  
**C3-7ED-N (2-5MHz) Convex, HL5-9ED-N Linear**  
**Like new very light use DEMO – ZONE A**
- 23 **SA-PICO Portable Digital Color Doppler (MFG: 2007) - NEW**  
Color Doppler, PW Doppler, M-Mode, OB/GYN/Vascular Pkg., SonoView, Tissue Harmonics, DICOM  
**C3-7ED-N Convex, L5-12ED-N (5-12MHz) Linear**  
**(1) Year Warranty – ZONE A**

### **SONOSITE**

- 24 **SonoSite TITAN (MFG: 2007), Remanufactured – Shared Service Application**  
Color Doppler, CW/PW Doppler, Angio, M-Mode, ECG, Mini Dock, THI  
**C15/4-2 Cardiac/Abd Micro Convex, L38/10-5 Vascular/Small Parts Linear**  
Operator's Manual and complete accessories  
**ZONE A**  
**6 month manufacturer warranty**
- 25 **SonoSite MicroMaxx (MFG: 2008), Remanufactured – Shared Service Application**  
Color Doppler, CW/PW Doppler, Angio, M-Mode, ECG, Mini Dock, THI, SonoRES  
**P17/5-1MHz Cardiac Sector, L38e/10-5 Vascular/Small Parts Linear**  
Operator's Manual and complete accessories – **like new in the BOX**  
**ZONE C**  
**5 years manufacturer warranty**

### **MINDRAY**

- 26 **MindRay DP-6600 (Portable) – MFG: 2007 – Reconditioned – DEMO UNIT**  
B/W system, USB Port for Image Storage on to Memory Stick, Dual Probe Port,  
10" Monitor, **35C50EA** – 3.0-6.0MHz Convex/50mmR  
**ZONE A**
- 27 **MindRay DP-3300 VET (Portable) – MFG: 2007 – Reconditioned – DEMO UNIT**  
B/W system, (2) USB Port for Image Storage on to Memory Stick  
10" Monitor, **65C15EAV** (5.0/6.5/8.0Mhz) Micro Convex  
**ZONE A**

### **CHISON**

- 28 **CHISON 8300 (Portable) – MFG: 2007 – NEW condition – Demo Unit**  
B/W System, USB image transfer, Harmonic Imaging, Dual Probe Port, 10" Monitor  
**C60613S (2.5-5.0MHz) Convex, L40617S (5.0-10.0MHz) Linear, C12616S (4.5-8.0MHz) Endo-vaginal**  
**ZONE A**
- 29 **CHISON 8300 (Portable) – MFG: 2007 – NEW condition – Demo Unit**  
B/W System, USB image transfer, Harmonic Imaging, Dual Probe Port, 10" Monitor  
**C60613S (2.5-5.0MHz) Convex, C12616S (4.5-8.0MHz) Endo-vaginal**  
**ZONE A**

## Stand-Alone Systems

### ACUSON (SIEMENS)

- 1 **SEQUOIA 512 LCD CVO – Shared Service – Rev. 12 (MFG: 2008)**  
**19" Flat panel LCD Monitor with articulating arm – Progressive Media & Display II Pkg.**  
**Cardiac/OB/Vascular Calc. Pkg. – with High Density Imaging Probe Option**  
Stress ECHO, Clarify VE, Coronary Flow Reserve, DICOM Store SCU/SR, DICOM Printing  
Dual Linear, Extended Trigger, GI Auto Doppler, GI Cadence, Hi Frame Rate DTI, High HR ECG, MIS/MIF  
Modality Worklist, Native TEQ, Spectral TEQ, Steered Spatial Compounding,  
Signature II, Driving Experience 7 & 8  
**4V1c** (1.5-4.0MHz) Adult Cardiac Sector, **9L4** Vascular/Small Parts Linear, **4C1** (1.5-4.0MHz) Convex  
**ZONE C**
- 2 **SEQUOIA 512 LCD CVO – Rev. 12. KPI 25137**  
**18" Flat panel LCD Monitor – Progressive Media & Display II Pkg.**  
**Cardiac/OB/Vascular Calc. Pkg.,** DICOM Bulk Store, DICOM Print, Modality Worklist,  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Stress Echo, ECG/AUX  
Calc Data to MO, Extended trigger, Duel Linear, Driving Experience 7 & 8, NTHI  
Native TEQ, High HR ECG, GI Auto Doppler, Steered Spatial Compounding  
Spectral TEQ, TEQ, Signature II Option  
**4V1c** (1.5-4.0MHz) Adult Cardiac Sector, **4C1** (1.5-4.0MHz) Convex, **8L5** (8-5Mhz) Linear  
**ZONE A**
- 3 **SEQUOIA 512 LCD – Rev. 12.**  
**18" Flat panel LCD Monitor – Progressive Media & Display II Pkg.**  
**OB/Vascular Calc. Pkg.,** DICOM Bulk Store, DICOM Print, Modality Worklist,  
Color Doppler, CW/PW Doppler, M-Mode, Angio, Calc Data to MO  
Duel Linear, Driving Experience 7 & 8, NTHI, Native TEQ  
GI Auto Doppler, Steered Spatial Compounding, Spectral TEQ, TEQ, Signature II Option  
**EV-8C4** (4-8Mhz) Transvaginal, **4C1** (1.5-4.0MHz) Convex, **8L5** (8-5Mhz) Linear  
**ZONE A**
- 4 **SEQUOIA 512 C/V/O (REV. 9.51) RE-MFG 2005 S/N: 62xxx KPI 24828**  
**Cardiac/OB/Vascular Calc. Pkg.,** Color Doppler, CW/PW Doppler, M-Mode,  
Angio, Stress Echo, ECG/AUX, Calc Data to Mo, **Dicom Bulk Store**, DTI,  
Driving Experience 7 & 8, Duel Linear, Extended Trigger, GI Auto Doppler,  
High Frame Rate DTI, **Modality Work list**, **NTHI**-Tissue Harmonics, Spectral TEQ, **TEQ**, Third Party Enabled  
**4V1c** (4-1Mhz) Cardiac, **6L3** (6-3Mhz) Vascular Linear, **4C1** (1.5-4.0MHz) Convex  
**ZONE A PENDING SALE**
- 5 **SEQUOIA 512 CVO Shared Service (Rev. 9.5) KPI 24057**  
NTHI -Tissue Harmonics, **DTI**, Driving Experience 7, Cadence  
**STRESS ECHO**, MOD, Color Doppler, CW/PW Doppler,  
M-Mode, Angio, ECG/AUX  
**3V2c** (2-3MHz) Cardiac Sector, **8L5** (8-5MHz) Linear, **4C1** (1.5-4.0MHz) Convex  
**ZONE B**
- 6 **SEQUOIA 512 CVO Shared Service Rev. 9.51 PAL/220V/50Hz (2005/2006 Software Version) KPI 22765**  
NTHI -Tissue Harmonics, 3<sup>rd</sup> Party Service Enabled, TEQ – Tissue Equalization  
DICOM Store/Print, DICOM Modality Worklist, MOD, STRESS ECHO, Signature II Option  
Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG/AUX, B-Color, Spatial Compounding  
**3V2c (2-3MHz)** Cardiac Sector, **8L5** (8-5MHz) Linear, **4C1** (1.5-4.0MHz) Convex  
**ZONE B**
- 7 **SEQUOIA 512 GI (Rev. 9.5) KPI 24065**  
OB/Vascular Calc. Pkg NTHI -Tissue Harmonics, MOD, Color Doppler, PW Doppler, Third Party Service  
M-Mode, Angio, **Dicom Store & Print**, Modality Worklist, Driving Experience 7 & 8, **TEQ**  
**EV-8C4** (4-8Mhz) Transvaginal, **8L5** (8-5MHz) Linear, **4C1** (1.5-4.0MHz) Convex, **4V1** (4-1Mhz) Vector  
**ZONE B**

- 8 **SEQUOIA 512 GI (Rev. 9.5) KPI 24827**  
OB/Vascular Calc. Pkg NTHI -Tissue Harmonics, MOD, Color Doppler, PW Doppler, M-Mode, Angio, **Dicom Store & Print**, Modality Worklist, Driving Experience 7, **TEQ EV-8C4** (4-8MHz) Transvaginal, **8L5** (8-5MHz) Linear, **4C1** (1.5-4.0MHz) Convex  
**ZONE B**
- 9 **SEQUOIA 512 GI (Rev. 9.5) KPI 24060**  
OB/Vascular Calc. Pkg NTHI -Tissue Harmonics, MOD, Color Doppler, PW Doppler, M-Mode, Angio, **Dicom Store**, GI Auto Doppler, Driving Experience 7 & 8, **Surface Rendering EV-8C4** (4-8MHz) Transvaginal, **8L5** (8-5MHz) Linear, **4C1** (1.5-4.0MHz) Convex  
**ZONE B**
- 10 **SEQUOIA C256 Cardiac/Vascular (Rev. 9.5)**  
STRESS ECHO, MOD, Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG/AUX, NTHI -Tissue Harmonics  
**3V2c** Adult Cardiac Sector, **8L5** (8-5MHz) Linear  
**ZONE B**  
"Special price"
- 11 **SEQUOIA C256 Cardiac/Vascular (Rev. 9.5)**  
NTHI -Tissue Harmonics, Modality Work list, Driving Experience 7 & 8, **Dicom Store / Print**, STRESS ECHO, MOD, Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG/AUX,  
**3V2c** Adult Cardiac Sector, **8L5** (8-5MHz) Linear  
**ZONE B**  
"Special price"
- 12 **SEQUOIA C256 Cardiac (Rev. 9.5)**  
NTHI -Tissue Harmonics, Modality Work list, Driving Experience 7 & 8, **Dicom Store**, STRESS ECHO, MOD, Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG/AUX, **3V2c** Adult Cardiac Sector  
**ZONE B**  
"Special price"
- 13 **SEQUOIA C256 Cardiac (Rev. 9.5)**  
NTHI -Tissue Harmonics, Cardiac Cadence 2, Driving Experience 7, STRESS ECHO, MOD, Color Doppler, CW/PW Doppler, **Dicom Store** M-Mode, Angio, ECG/AUX, **3V2c** Adult Cardiac Sector  
**ZONE B**  
"Special price"

## **ATL (PHILIPS)**

- 14 **UM-400C Color Doppler (Rev. 2.05) – OB/Vascular Calc. (MFG: 2000-2001) – upgraded 2008**  
Color Doppler, PW Doppler, M-Mode, Power Doppler, Surface 3D  
**CLA3.5/40mm** Convex, **EC6.5** Endo-cavity  
**ZONE A**
- 15 **UM-400C Color Doppler (Rev. 2.05) – OB/Vascular Calc. (MFG: 2001-2002) - upgraded 2008**  
Color Doppler, PW Doppler, M-Mode, Power Doppler, Surface 3D, Tissue Harmonics  
**CLA4.0/40mm** Convex, **EC6.5** Endo-cavity  
**ZONE A -- NEW PROBES**
- 16 **UM-400C Color Doppler (Rev. 2.05) – OB/Vascular Calc. (MFG: 2002) – upgraded 2008**  
Color Doppler, PW Doppler, M-Mode, Power Doppler, Surface 3D, Tissue Harmonics  
**C4.0** Convex, **EC6.5** Endo-cavity, **HL7.5** (5-9MHz) Linear  
**ZONE A -- NEW PROBES**
- 17 **ATL HDI-4000 Live 3D MFG: 2002**  
Color Doppler, PW Doppler, M-Mode, Angio, CD Writer, MOD  
Live 3D/4D, Volume 3D, Surface 3D, SonoView, SonoAtlas, CAFÉ, Tissue Harmonics, Pulse Inversion  
**3D7-4** (Real-time 4D) Convex, **3D8-5v** Volume 3D Endo-vaginal, **C5-2** Curved Linear  
**ZONE C**

- 18 **PHILIP HDI-4000 Live 3D**  
Color Doppler, PW Doppler, M-Mode, Angio, CD Writer, MOD  
Live 3D/4D, Volume 3D, Surface 3D, SonoView, SonoAtlas, CAFÉ, Tissue Harmonics, Pulse Inversion  
**3D7-4** (Real-time 4D) Convex, **3D8-5v** Volume 3D Endo-vaginal  
**SOLD**
- 19 **HDI-5000 OB/Vascular (Rev. 127.xx)**  
**Color Doppler**, Power Angio, PW Doppler, M-mode, **Tissue Harmonic Application**,  
Net Link, Disk Link, DICOM, TSI, Hi-Q,  
**C5-2** Curved Linear, **L7-4** Vascular Linear  
**ZONE A**
- 20 **HDI-5000 SonoCT OB/Vascular (Rev. 190.xx)**  
**2<sup>nd</sup> Gen. SonoCT (Compound Imaging), Color Doppler**, Power Angio, PW Doppler  
M-mode, **Tissue Harmonic Application**,  
Net Link, Disk Link, DICOM, TSI, Hi-Q,  
**C5-2** Curved Linear, **CL10-5 Compact** Linear (Hockey Stick)  
**ZONE A**
- 21 **HDI-5000 SonoCT OB/Vascular (Rev. 190.xx)**  
**2<sup>nd</sup> Gen. SonoCT (Compound Imaging), Color Doppler**, Power Angio, PW Doppler  
M-mode, **Tissue Harmonic Application**,  
Net Link, Disk Link, DICOM, TSI, Hi-Q,  
**C5-2** Curved Linear, **L7-4** Vascular Linear  
**ZONE A**
- 22 **HDI-5000 SonoCT OB/Vascular (Rev. 190.xx)**  
**2<sup>nd</sup> Gen. SonoCT (Compound Imaging), Color Doppler**, Power Angio, PW Doppler, M-mode, 3D Angio  
**Tissue Harmonic Application**, Net Link, Disk Link, DICOM, TSI, Hi-Q,  
**C5-2** Curved Linear, **L7-4** Vascular Linear  
**ZONE A**
- 23 **HDI-5000 SonoCT OB/Vascular (Rev. 190.xx)**  
**2<sup>nd</sup> Generation SonoCT, Color Doppler**, Power Angio, CW/PW Doppler, M-mode, 3D Angio  
**Tissue Harmonic Application**, Net Link, Disk Link, DICOM, TSI, Hi-Q,  
**C5-2** Curved Linear, **C8-4v (4-8MHz)** Endo-vaginal  
**ZONE A**
- 24 **HDI-5000 SonoCT CVO (Rev. 190.xx) – Shared Service**  
**2nd Gen. Sono-CT**, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
**Tissue Harmonic Application**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q  
**C5-2** Curved Linear, **L7-4** Vascular Linear, **P4-2** Phased Array Cardiac  
**ZONE A**
- 25 **HDI-5000 SonoCT CVO (Rev. 190.xx) – Shared Service**  
**2nd Gen. Sono-CT**, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
**Tissue Harmonic Application**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q  
**C5-2** Curved Linear, **L12-5/38mm** Vascular, **C8-4v** Endo-vaginal  
**ZONE A**
- 26 **HDI-5000 SonoCT CVO (Rev. 198.xx) – Shared Service**  
**2nd Gen. Sono-CT**, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
**Tissue Harmonic Application**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q  
**C5-2** Curved Linear, **L12-5/38mm** Vascular, **C8-4v** Endo-vaginal, **L7-4** Vascular Linear  
**ZONE A**
- 27 **HDI-5000 SonoCT CVO (Rev. 210.xx) – Shared Service**  
**3rd Gen. Sono-CT**, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
Panoramic Imaging, **Tissue Harmonic**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q  
DICOM Worklist, **C5-2** Curved Linear, **L7-4** Linear, **P4-2** Cardiac Sector  
**ZONE B**

- 28 **HDI-5000 SonoCT XRES CVO (Rev. 210.xx) – Shared Service**  
3rd Gen. Sono-CT, XRES, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
Panoramic Imaging, **Tissue Harmonic**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q  
DICOM Worklist, Advanced 3D, **C5-2** Curved Linear, **C8-4v** Endo-vaginal  
**ZONE B**
- 29 **HDI-5000 SonoCT XRES CVO (Rev. 210.xx) – Shared Service**  
3rd Gen. Sono-CT, XRES, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
Panoramic Imaging, **Tissue Harmonic**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q  
DICOM Worklist, Advanced 3D,  
**C5-2** Curved Linear, **L7-4** Vascular Linear , **P4-2** Cardiac Sector  
**ZONE B**
- 30 **PHILIPS 5000 SonoCT XRES CVO (Rev. 210.xx) – Shared Service (MFG: 2003)**  
3rd Gen. Sono-CT, XRES, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
Panoramic Imaging, **Tissue Harmonic**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q,  
**Panasonic SVHS VCR**, **Sony UP-895MD** B/W Printer, DICOM Worklist, Advanced 3D, HDI Speed-up  
**C5-2** Curved Linear, **L7-4** Vascular Linear, **L12-5/38mm** Vascular, **C8-4v** Endo-vaginal, **P4-2** Cardiac Sector  
**ZONE A**
- 31 **PHILIPS 5000 SonoCT XRES CVO (Rev. 210.xx) – Shared Service (MFG: 2003)**  
3rd Gen. Sono-CT, XRES, Power Angio, Color Doppler, CW/PW Doppler, M-mode, ECG/AUX  
Panoramic Imaging, **Tissue Harmonic**, Net Link, Disk Link, Scaled DICOM, TSI, Hi-Q,  
**Panasonic SVHS VCR**, **Sony UP-895MD** B/W Printer, DICOM Worklist, Advanced 3D, HDI Speed-up  
**C5-2** Curved Linear, **L12-5/38mm** Vascular, **C8-4v** Endo-vaginal, **P4-2** Cardiac Sector  
**ZONE A -- NEW PROBES**

## **Philips/Agilent/HP**

- 32 **HP Sonos 1000 CVO REV 0.4 (MFG 1998)**  
Shared Service Application: Cardiac/OB/Vascular, B-Color  
Color Doppler, CW/PW Doppler, M-Mode, Angio, ECG/AUX, Panasonic AG-7350 SVHS VCR  
**21202A** (2-2.5Mhz) Cardiac sector  
**ZONE B**
- 33 **Sonos 5500 STRESS, Color Doppler, Rev. B.2.3, Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 512MB Memory, MCE, HI MCE, Stress-Echo, Contrast  
Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR  
**S3 (21311A)** Adult Cardiac Sector  
**ZONE B**
- 34 **Sonos 5500 STRESS, Color Doppler, Rev. B.2.3, Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 512MB Memory, MCE, HI MCE, Stress-Echo, Contrast  
Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR  
**S3 (21311A)** Adult Cardiac Sector, **L5040** (21355B) Vascular Linear  
**ZONE B**

### **Optional Transducers for Sonos 5500/ 7500: (not for individual sales)**

- 21321A (C3540) Ultra-band Convex: recon/new
- 21336A (E6509) Endo-cavity: recon/new
- 21355B (L5040) 3.0-7.0MHz Trap Linear: (Rev. A & B)
- 21390A (15L6) Hockey Stick Linear: recon/new
- 21356A (11-3L) linear: recon/new
- 21358B (L7540) Linear:
- 21350A (S8) Pediatric Sector
- 21380A (S12) Pediatric Sector
- 21221B (1.9CW) 21221A (1.9CW)

#### **TEE PROBES:**

- 21369A (OMNI II) TEE
- 21378A (OMNI 3) TEE

#### **OPTIONS: (Must order with the system)**

- 1) IDI-DICOM (for Rev.B.0 +)
- 2) SonosPC DICOM with CD Writer (for Rev. D.1 +)

- 35 **Sonos 5500 STRESS, Color Doppler, Rev. D.2, Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 2GB Memory, MCE, HI MCE, Stress-Echo, Contrast Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR **S3 (21311A) Adult Cardiac Sector**  
**ZONE B**
- 36 **Sonos 5500 STRESS, Color Doppler, Rev. D.2, Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 2GB Memory, MCE, HI MCE, Stress-Echo, Contrast Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR **S3 (21311A) Adult Cardiac Sector, 21369A (OMNI II) TEE**  
**ZONE B**
- 37 **Sonos 5500 STRESS, Color Doppler, Rev. D.2, Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 2GB Memory, MCE, HI MCE, Stress-Echo, Contrast Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR **S3 (21311A) Adult Cardiac Sector, 15L6 (21390A) 6-15MHz Hockey Stick Linear**  
**SPECIAL PKG. Price – ZONE B**
- 38 **Sonos 5500 STRESS, Color Doppler, Rev. D.2, Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 2GB Memory, MCE, HI MCE, Stress-Echo, Contrast Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR **S3 (21311A) Adult Cardiac Sector, 11-3L (21356A) Small Parts/Vascular Linear (21356A Linear – NEW) – ZONE B**
- 39 **Sonos 5500 STRESS, Color Doppler, Rev. D.2, Stress-Echo, Contrast – Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 2GB Memory, MCE, HI MCE, **SONOS PC DICOM** Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR **S3 (21311A) Adult Cardiac Sector, 11-3L (21356A) Small Parts/Vascular Linear (21356A Linear – NEW) – ZONE B**
- 40 **Sonos 5500 STRESS, Color Doppler, Rev. D.2, Stress-Echo, Contrast – Shared Service (CVO)**  
Fusion Harmonics, Study Manager, DSR, 2GB Memory, MCE, HI MCE, **SONOS PC DICOM** Color Doppler, Doppler Tissue, B-Color, CW/PW Doppler, ECG/AUX, Panasonic SVHS VCR **S3 (21311A) Adult Cardiac Sector, 11-3L (21356A) Small Parts/Vascular Linear, C3540 (2.5-6.0MHz) Convex (21356A Linear & C3540 (21321A) – NEW probes) – ZONE B**
- 41 **PHILIPS SONOS 7500 Rev. D.2 (MFG: 2003) KPI 24779**  
**SonosPC DICOM, Fusion Harmonics, Stress Echo, Contrast, Study Manager, 2GB Memory, MCE, HI MCE** Color Doppler, CW/PW Doppler, M-Mode, Angio, B-Color, Tissue Doppler, ECG/AUX **S3 (21311A) Adult Cardiac Sector, Panasonic SVHS VCR**  
**ZONE A**
- 42 **PHILIPS SONOS 7500 Live 3D Rev. D.2 (MFG: 2004) KPI 24778**  
**Fusion Harmonics, Stress Echo, Contrast Imaging**  
**LIVE 3D, Real-time Bi-plane Imaging SonosPC DICOM, Study Manager, 2GB Memory, MCE, HI MCE, Color Doppler, CW/PW Doppler, M-Mode, Angio, B-Color, Tissue Doppler, ECG/AUX**  
**S3 (21311A) Adult Cardiac Sector, 11-3L (21356A) Small Parts/Vascular Linear (NEW), Panasonic SVHS VCR**  
**ZONE A**
- 43 **PHILIPS SONOS 7500 Live 3D Rev. D.2 (MFG: 2003) KPI 24777**  
**LIVE 3D, Realtime Bi-plane Imaging, Fusion Harmonics, Stress Echo, Contrast Imaging**  
**SonosPC DICOM, Study Manager, 2GB Memory, CD Writer, MCE, HI MCE** Color Doppler, CW/PW Doppler, M-Mode, Angio, B-Color, Tissue Doppler, ECG/AUX **X4 (21315A) Live 3D Cardiac Sector, S3 (21311A) Adult Cardiac Sector, Panasonic SVHS VCR**  
**ZONE A**

## **PHILIPS**

- 44 **PHILIPS HD-3 OB/Vascular (MFG: 2008) Like new in the box**  
**Color Doppler, M-Mode, Angio, Tissue Harmonics, LCD Flat Panel Monitor, DICOM** Surface 3D, SonoView  
**C3-7 Curved “New Probe”**  
**ZONE A**

- 45 **PHILIPS HD-3 OB/Vascular (MFG: 2008) Like new in the box**  
Color Doppler, M-Mode, Angio, Tissue Harmonics, LCD Flat Panel Monitor, DICOM  
Surface 3D, SonoView  
C3-7 Curved, EC4-9 Endo-cavity "New Probes"  
**ZONE A**
- 46 **PHILIPS HD-3 OB/Vascular (MFG: 2008) Like new in the box**  
Color Doppler, M-Mode, Angio, Tissue Harmonics, LCD Flat Panel Monitor, DICOM  
Surface 3D, SonoView  
C3-7 Curved, L5-9 Linear, EC4-9 Endo-cavity "New Probes"  
**ZONE A**
- 47 **PHILIPS ENVISOR HD OB/Vascular (Rev. C.1.3) – MFG: 2003 -- upgraded in 10/2008**  
OB/Vascular Pkg. Color Doppler, CW/PW Doppler, M-Mode, CD Writer  
DICOM Networking, Performance Pkg., i-Scan  
C5-2 (21426A) Curved Linear, L12-5 Small Parts/Vascular Linear  
**ZONE A**
- 48 **PHILIPS ENVISOR HD CVO (Rev. C.1.3) – MFG: 2003 -- upgraded in 10/2008**  
Cardiac/OB/Vascular Pkg. Color Doppler, CW/PW Doppler, M-Mode, CD Writer  
DICOM Networking, Performance Pkg, i-Scan  
C5-2 (21426A) Curved Linear, L1038 (21376A) Small Parts/Vascular Linear, C8-4v (21337A) Endovaginal  
PA4-2 (21422A) Cardiac Sector  
**ZONE A**
- 49 **PHILIPS ENVISOR HD CVO (Rev. C.1.3) – MFG: 2003 - upgraded in 10/2008**  
Cardiac/OB/Vascular Pkg. Color Doppler, CW/PW Doppler, M-Mode, CD Writer  
DICOM Networking, Performance Pkg., i-Scan  
PA4-2 (2-4MHz) Cardiac Sector, 21369A (OMNI II) TEE, L7535 (21359A) Linear  
**ZONE A**
- 50 **PHILIPS HD-11XE OB/Vascular (MFG: 2006-2007)**  
17" LCD Monitor, Vascular/OB/GYN/General Imaging Pkg.  
Color Doppler, Steered CW, PW Doppler, M-Mode, Angio, CD Writer  
SonoCT, XRES, i-SCAN, DICOM Structured Reporting, DICOM Networking  
C6-3 Curved Linear, L8-4 Vascular Linear  
**ZONE A**
- 51 **PHILIPS HD-11XE OB/Vascular (MFG: 2007) – Like New in the BOX**  
17" LCD Monitor, Vascular/OB/GYN/General Imaging Pkg.  
Color Doppler, Steered CW, PW Doppler, M-Mode, Angio, CD Writer  
SonoCT, XRES, i-SCAN, DICOM Structured Reporting, DICOM Networking  
C5-2 Curved Linear, C8-4v Endo-vaginal, L12-3 Linear  
**ZONE A**
- 52 **PHILIPS HD-11XE CVO (MFG: 05/2007) – Like New in the BOX**  
17" LCD Monitor, Shared Service, Cardiac/Vascular/OB/GYN/General Imaging Pkg.  
Color Doppler, Steered CW, PW Doppler, M-Mode, Angio, ECG/AUX, CD Writer  
SonoCT, XRES, i-SCAN, DICOM Structured Reporting, DICOM Networking  
S3-1 (Cardiac) Sector, C6-3 Curved Linear, C8-4v Endo-vaginal, L12-3 Linear  
**ZONE A**
- 53 **PHILIPS HD-11XE CVO (MFG: 05/2007) – Like New in the BOX**  
17" LCD Monitor, Shared Service, 3D/4D Option, Cardiac/Vascular/OB/GYN/General Imaging Pkg.  
Color Doppler, Steered CW, PW Doppler, M-Mode, Angio, CD Writer  
SonoCT, XRES, i-SCAN, DICOM Structured Reporting, DICOM Networking  
C5-2 Curved Linear, C8-4v Endo-vaginal, 3D8-4 3D Convex  
**ZONE A**
- 54 **PHILIPS HD-11XE C/V– STRESS ECHO (MFG: 04/2007) – Like New in the BOX**  
17" LCD Monitor, Shared Service, Cardiac/Vascular/General Imaging Pkg.  
Color Doppler, Steered CW, PW Doppler, M-Mode, Angio, ECG/AUX, CD Writer, Stress Echo  
SonoCT, XRES, i-SCAN, DICOM Structured Reporting, DICOM Networking  
S3-1 (Cardiac) Sector, C6-3 Curved Linear, L12-3 Linear  
**ZONE A**

**KPI Ultrasound, Inc.**  
**887 Marlborough Ave**  
**Riverside, CA 92507**

**Tel: (951) 367-0872**  
**Fax: (951) 367-0891**  
**[www.kpiultrasound.com](http://www.kpiultrasound.com)**

**iU-22 / iE-33 Available Transducers in Stock:**

C5-2 (5-2MHz) Curved Linear – **SOLD**  
C9-4 (9-4MHz) Curved Linear -  
C9-5ec (9-4MHz) Endo-cavity -  
C8-4v (8-4MHz) Endo-vaginal - **SOLD**  
L17-5 (17-5MHz) Hi-Frequency Linear – new  
L15-7io (15-7MHz) Hockey Stick Linear – new  
L12-5 (12-5MHz) Linear – **SOLD**  
11-3L (11-3MHz) Linear -- (for iE-33 only) -  
L9-3 (9-3MHz) Linear - new  
L8-4 (8-4MHz) Linear -  
S5-1 (5-1MHz) PureWave Cardiac Sector – **SOLD**  
S3-1 (3-1MHz) Cardiac Sector (for iU-22)-  
S4-1 (4-1MHz) Deep Abdominal Imaging -  
S8-3 (8-3MHz) Pediatric Cardiac Sector -  
X3-1 (1-3MHz) Live 3D cardiac Sector – demo  
X7-2 (7-2MHz) Live 3D Pediatric cardiac sector - demo  
3D6-2 (6-2MHz) 3D/4D Volume Convex (for iU-22/HD-11XE) -  
3D8-4 (8-4MHz) 3D/4D Volume Convex (for HD-11/XE) -

**TEE Probe**

S7-3t (7-3MHz) Mini Multi TEE - (Demo)

**GE (General Electric)**

- 55 **LOGIQ 200 PRO (MFG: 2004)** B/W System  
ATO – Auto Tissue Optimization, Comprehensive OB Calc.  
3CB (2-5MHz) Convex, MTZ (4-9MHz) Endo-cavity  
**ZONE A**
  
- 56 **LOGIQ 200 PRO (MFG: 2002)** B/W System  
ATO – Auto Tissue Optimization, Comprehensive OB Calc.  
3CB (2-5MHz) Convex, MTZ (4-9MHz) Endo-cavity  
**ZONE A “NEW TRANSDUCERS”**
  
- 57 **LOGIQ 400 MR3 REV4.01 (MFG: 1997)** KPI 24683  
Color Doppler, PW Doppler, M-Mode, Angio, Auto Doppler, THI  
C358 (2-5MHz) Convex, E721 (4-9MHz) Endo-cavity  
**ZONE B PENDING SALE**
  
- 58 **LOGIQ 3 PRO BT05 (MFG: 2003) Rev. 4.1.2**  
Cardiac/OB/GYN/Small Part/Urology, Multi-Frequency, ATO (Auto Tissue Optimization), THI,  
Steered CW, Anatomic M-Mode,  
3.5C (2.0-5.0MHz) Convex, E8C Endovaginal  
**ZONE B**
  
- 59 **LOGIQ 3 PRO BT05 (MFG: 2005) Rev. 4.0.0**  
OB/GYN/Small Part/Urology, Multi-Frequency, ATO (Auto Tissue Optimization), THI, DICOM  
3.5C (2.0-5.0MHz) Convex, E8C Endovaginal  
**ZONE B**
  
- 60 **LOGIQ 3 PRO BT05 (MFG: 8/2005) Rev. 4.0.0**  
OB/GYN/Small Part/Urology, Multi-Frequency, ATO (Auto Tissue Optimization)  
Steered CW, Anatomical M-Mode (AMM), Easy 3D, DICOM  
3S (1.5-3.5MHz) Cardiac Sector, 3.5C (2.0-5.0MHz) Convex, 10LB (5-10MHz) Small Parts/Vascular Linear  
**ZONE C**
  
- 61 **LOGIQ 5 PRO BT04 (MFG: 2004) KPI 24872**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO, Virtual Convex, THI, CD Writer  
3.5C (2-5MHz) Convex, E8C (4-10MHz) Endo-cavity  
**ZONE B**

- 62 **LOGIQ 5 PRO BT05 (MFG: 2005) KPI 24871**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO, Virtual Convex, THI, CD Writer , LogiqView, Anatomic M-Mode  
**3.5C** (2-5MHz) Convex, **E8C (4-10MHz)** Endo-cavity  
**ZONE B**
- 63 **LOGIQ 5 Expert BT05 OB/VASC (MFG: 2005) KPI 24833**  
Color Doppler, PW Doppler, M-Mode, Angio, Dicom, ATO, Virtual Convex, THI, CD Writer  
**3.5C** (2-5MHz) Convex, **10L** (5-10MHz) Small Parts/Vascular Linear  
**ZONE B - NEW PROBES**
- 64 **LOGIQ 5 Expert BT05 CVO (MFG: 2007) KPI 24869**  
Color Doppler, CW/ PW Doppler, M-Mode, Angio, ATO, Virtual Convex, THI, CD Writer, ECG  
**3.5C** (2-5MHz) Convex, **3S** (1.5-3.5MHz) Cardiac Sector, **10L** (5-10MHz) Small Parts/Vascular Linear  
**ZONE B - NEW PROBES**
- 65 **LOGIQ P5 (MFG: 12/2006) - 15" TFT LCD Monitor KPI 24870**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO, Logiq VUE  
Virtual Convex, THI, CD Writer, OB/Vasc/Card Calc.  
**3.5C** (2-5MHz) Convex, **E8C (4-10MHz)** Endo-cavity  
**ZONE A**
- 66 **LOGIQ P5 (MFG: 12/2006) - 15" TFT LCD Monitor / 3D-4D Package**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO, Logiq VUE  
Virtual Convex, THI, CD Writer, OB/Vasc/Card Calc.  
**3.5C** (2-5MHz) Convex, **E8C (4-10MHz)** Endo-cavity, **4D3C-L** (2.0-5.0Mhz) 4D Convex  
**ZONE A**
- 67 **LOGIQ S6 BTO6 CVO Shared Service (MFG:2007) W/ 4D Imaging KPI 24897**  
**LCD Monitor**, Color Doppler, CW/PW Doppler, M-Mode, Angio, ATO,  
Virtual Convex, THI, CD Writer, **4D Imaging**, Cross Beam, ECG  
**4C (2-5MHz)** Convex, **E8C** (4-9MHz) End-cavity, **3S** (1.5-3.5MHz) Cardiac Sector  
**12L** (4.5-13Mhz) Vascular/Small Parts  
**ZONE A**
- 68 **LOGIQ 9 BT07 LCD Monitor with Real-time 4D KPI 24934**  
Color Doppler, PW Doppler, Power Doppler, M-Mode, THI, ATO,  
RT4D, Speckle Reduction (SRI), Crossbeam (CRI), DICOM, LOGIQVIEW, B-Flow, PFD,  
**4D3C (2.0-5.0MHz)** 4D convex, **E8C** (4-9MHz) End-cavity, **4C (2-5MHz)** Convex, **M12L** (4.5-13Mhz) Small Parts  
**ZONE A PENDING SALE**
- 69 **Vivid 7 Expert Dimension BT08 LCD Monitor (MFG 2007)**  
**20" Flat Panel Display** Color Doppler, CW/PW Doppler, M-Mode, Angio, Dicom Media,  
Anat. M-Mode, TVI Tissue tracking, Q Scan w/ Strain, Q Analysis, Blood Flow Image, Echo Pac, Stress Echo, Mpeg  
Vue, Echo Dicom, Modality Worklist, Advanced Q Scan w/ TSI, Multi-Dimensional Imaging w/ Bi-plan  
ECG/AUX, DVD R/W, MOD 9.1  
**M4S** Cardiac Sector, **9L** (3-7MHz) Vascular/Small parts Linear  
**ZONE A PENDING SALE**
- 70 **Vivid 7 Expert BT03 (MFG 2003)**  
Color Doppler, CW/PW Doppler, M-Mode, Angio, **Stress echo**, Tissue tracking, Archive, **Dicom**  
Plug and Scan, Q Analysis, ECG/AUX, DVD R/W, MOD  
**M3S** Cardiac Sector, **7L** (3-7MHz) Vascular/Small parts Linear  
**ZONE B**
- 71 **VOLUSON 730 PRO BT02 (MFG: 2003) - K23464**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, Tissue Harmonics, CD Writer, MOD  
**Real-time 4D, Volume 3D, Surface 3D, Multiplane 3D, DICOM, VOCAL**  
**RAB4-8P** (4-8MHz) 4D Convex  
**ZONE A**
- 72 **VOLUSON 730 PRO BT02 (MFG: 2003) K24641**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, Tissue Harmonics, CD Writer, MOD  
**Real-time 4D, Volume 3D, Surface 3D, Multiplane 3D**  
**RAB2-5** (2-5MHz) 4D Convex, **AB2-7** (2-7MHz) 2D Convex, **SP6-12** Vascular/Small parts Linear  
**ZONE A**

- 73 **VOLUSON 730 PRO-V BT04**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, CD Writer  
**Real-time 4D, Volume 3D, Surface 3D, Multiplane 3D**  
**RAB4-8P** (2-5MHz) 4D Convex, **AB2-7** (2-7MHz) 2D Convex  
**ZONE A PENDING SALE**
- 74 **VOLUSON 730 PRO-V BT04 (MFG: 2005)**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, CD Writer  
**Real-time 4D, Volume 3D, Surface 3D, Multiplane 3D**  
**RAB2-5L** (2-5MHz) 4D Convex, **RIC5-9H** (5-9MHz) 4D Vaginal  
**ZONE C**
- 75 **VOLUSON 730 Expert BT02 (REV 2.2.1)**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, CD Writer  
4D Imaging, **DICOM, VOCAL, STIC**, Harmonics  
**RAB4-8P** (2-5MHz) 4D Convex, **RIC5-9** (5-9MHz) 4D Vaginal  
**ZONE B**
- 76 **VOLUSON 730 Expert BT02 (REV 2.2.1) KPI 24894**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, CD Writer  
4D Imaging, **DICOM, CRI, VCI**  
**RAB4-8P** (2-5MHz) 4D Convex, **RIC5-9** (5-9MHz) 4D Vaginal  
**ZONE A PENDING SALE**
- 77 **VOLUSON 730 Expert BT05 (MFG: 02/2006)**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, Tissue Harmonics, DVD Writer  
**Real-time 4D, Volume 3D, Surface 3D, DICOM, SonoView II, (CRI Crossbeam)**  
**SRI II (Speckle Reduction), TUI (Tomographic Ultrasound Imaging) – Multi-Slice**  
**B-Flow, STIC**  
**RAB4-8L** (4-8Mhz) 4D Convex, **RIC5-9H** (5-9Mhz) 4D Vaginal, **4C-A** (2-5Mhz) Convex  
**ZONE A PENDING SALE**
- 78 **VOLUSON 730 Expert BT08 (MFG: 05/2006)**  
Color Doppler, PW Doppler, M-Mode, Angio, ATO- Auto Tissue, Tissue Harmonics, DVD Writer  
**Real-time 4D, Volume 3D, Surface 3D, DICOM, SonoView II, (CRI Crossbeam)**  
**SRI II (Speckle Reduction), TUI (Tomographic Ultrasound Imaging) – Multi-Slice**  
**B-Flow, STIC**  
**RAB4-8L** (4-8Mhz) 4D Convex, **RIC5-9H** (5-9Mhz) 4D Vaginal, **4C-A** (2-5Mhz) Convex  
**ZONE B (NEW PROBES)**

## **MEDISON**

- 79 **SA-PICO Portable Digital Color Doppler (MFG: 2007) – DEMO UNIT**  
**Color Doppler**, PW Doppler, M-Mode, OB/GYN/Vascular Pkg., SonoView, Surface 3D, DICOM  
**HC2-5ED-N** (2-5MHz) Convex, **EC4-9ED-N** (4-9MHz) Endocavity  
**ZONE A**
- 80 **SA-PICO Portable Digital Color Doppler (MFG: 2007) --- NEW with latest version 2008**  
**Color Doppler**, PW Doppler, M-Mode, OB/GYN/Vascular Pkg., SonoView, Surface 3D, DICOM  
**C3-7ED-N** Convex, **L5-12ED-N (5-12MHz) Small Parts & Vascular Linear**  
**ZONE A**
- 81 **MEDISON SA-9900MT PRIME with Live 3D (upgraded 2008) – Like new condition – Rev. 2.08**  
Real-time Live 3D, Multi-planar 3D, Surface 3D – Premium 3D Pkg.  
Color Doppler, PW Doppler, Power Doppler, M-Mode  
Tissue Harmonics, Pulse Inversion Harmonic, DICOM, CD Writer, MOD  
**3D4-7EK Real-time Live 3D Convex, C3-7IM (2D Convex), EC4-9IS (2D Endocavity)**  
**ZONE C W/ NEW PROBES**
- 82 **SA-8000LV Prime – Live 3D (RE-MFG: 2008) – “Like New”**  
Color Doppler, PW Doppler, Angio, M-Mode, Tissue Harmonics, Pulse Inversion Harmonics  
Live 3D/Multi-Plane 3D/Surface 3D, SonoView, DICOM, CAFÉ, CD Writer  
**3D4-7EK (4-7MHz) 3D Convex, C3-7ED-N (3-7MHz) 2D Convex, EC4-9ED-N (4-9MHz) 2D Endo-Vag.**  
**ZONE C W/ NEW PROBES**

- 83 **SA-8000LV Prime – Live 3D (Reconditioned 2008)**  
Color Doppler, PW Doppler, Angio, M-Mode, Tissue Harmonics, Pulse Inversion Harmonics  
Live 3D/Multi-Plane 3D/Surface 3D, SonoView, DICOM, CAFÉ, CD Writer  
3D4-7EK (4-7MHz) 3D Convex  
**ZONE C**
- 84 **MEDISON ACCUVIX XQ (MFG: 2004) (Reconditioned 2008) – Rev. 1.07**  
Real-time Live 3D, Multi-planar 3D, Surface 3D – Premium 3D Pkg.  
Color Doppler, PW Doppler, Power Doppler, M-Mode  
Tissue Harmonics, Pulse Inversion Harmonic  
3DXI (Multi-Slice Tomographic Image), DICOM  
3D7-4EK Real-time Live 3D Convex, C3-7IM 2D Convex, L5-12IM 2D Linear  
**ZONE A**
- 85 **MEDISON ACCUVIX XQ (MFG: 2007) – Demo Unit – Rev. 1.07 - Like NEW**  
Real-time Live 3D, Multi-planar 3D, Surface 3D – Premium 3D Pkg.  
Color Doppler, PW Doppler, Power Doppler, M-Mode  
Tissue Harmonics, Pulse Inversion Harmonic  
3DXI (Multi-Slice Tomographic Image), DYNAMIC MR, DICOM  
3D4-8ET Real-time Live 3D Convex (NEW), C3-7IM (2D Convex), EC4-9IS (2D Endocavity)  
**ZONE A**

## **TOSHIBA**

- 86 **SSA-390A (PowerVision 8000, Rev.3.0), MFG: 2000 -2001 - Shared Service C/V/O**  
512 Channel, Color Doppler, CW/PW Doppler, Angio, M-Mode, Tissue Harmonics  
ECG/AUX, Tissue Harmonics, Tissue Doppler, DICOM (PowerView)  
PLN-503AT (4.0-8.0MHz) Vascular Linear, PSN-25AT  
**ZONE C**
- 87 **SSA-390A (PowerVision 8000, Rev.3.0), MFG: 2000 -2001 - Shared Service C/V/O**  
512 Channel, Color Doppler, CW/PW Doppler, Angio, M-Mode, Tissue Harmonics  
ECG/AUX, Tissue Harmonics, Tissue Doppler, DICOM (PowerView)  
PSN-25AT, PLN-703AT (4.0-8.0MHz) Vascular Linear,  
**ZONE C**
- 88 **SSA-390A (PowerVision 8000, Rev.3.0), MFG: 2000 -2001 - Shared Service C/V/O**  
512 Channel, Color Doppler, CW/PW Doppler, Angio, M-Mode, Tissue Harmonics  
ECG/AUX, Tissue Harmonics, Tissue Doppler, DICOM (PowerView)  
PVN375AT Convex, PLN-805AT (6.0-12.0MHz) Linear for Breast and Small Parts  
**ZONE C**
- 89 **SSA-550A (NEMIO 20) Color Doppler – MFG: 2001 (Reconditioned 2008)**  
OB/GYN/Vascular Application, Color Doppler, PW Doppler, DICOM, THI, ETUDE  
PVM-375AT (3.0-6.0MHz) Convex, PLM-703AT Vascular Linear  
**ZONE B**
- 90 **SSA660A (Xario Rev. 5.0) MFG: 2005 – OB/Vasc**  
Radiology Pkg. Color Doppler, PW Doppler, Power Doppler, M-Mode, Tissue Harmonics Imaging, ApliPure,  
I-Assist, Basic Dicom, CD Writer, Migration Pkg, Trapezoid Imaging, Quick Scan, B/W Sony UP-895MD Printer  
PVT-375BT (3-6MHz) Convex (NEW), PLT-704AT (6-11MHz) Linear  
**ZONE A**

## **SIEMENS**

- 91 **Siemens Antares (MFG: 2007) REV. 5.0**  
LCD Flat Panel with articulating arm, OB/GYN/Vascular Pkg. Color Doppler, PW Doppler  
M-Mode, THI, DIMAQ, Image Data Management Sieclear, Summit, Syngo General License, Film VFS,  
TEQ, GET Worklist, Film Paper Printer, Clarify VE, Boost, Doppler TEQ, MPPS Support,  
Advanced Sieclear SC, Security Privacy Basic  
Ch4-1(4-1Mhz) Convex, EC9-4 (9-4MHZ) OB/GYN-Urology, VF13-5 (13-5Mhz) Linear, VFX9-4 (9-4Mhz) Linear  
**ZONE A PENDING SALE**

- 92 **Siemens X300 LCD (MFG: 2007-2008) New in the Original Box**  
15" Flat Panel Display, Cardiac/OB/GYN/Vascular Pkg. Color Doppler  
CW/PW Doppler, M-Mode, THI, DIMAQ, Dicom, Image Data Management, CD-Writer, B/W printer  
**CH5-2** (2-5Mhz)Convex, **VF10-5** (5-10Mhz)Linear, **P4-2** (2-4Mhz) Cardiac , **2.0Mhz** PEDOF  
**ZONE C**
- 93 **Siemens G40 Color Doppler (MFG: 2008) Rev. 2.0 - NEW SYSTEM**  
OB/GYN/Vascular Pkg. Color Doppler, PW Doppler, M-Mode, THI, DIMAQ, Image Data Management  
CD-Writer, Mitsubishi P93 B/W printer, DICOM, **CH5-2** Convex, **EC9-4** Endo-cavity  
**w/ (1) Year Warranty -- 220/110V -- ZONE C**
- 94 **Siemens G40 Color Doppler (MFG: 2008) Rev. 2.0 - NEW SYSTEM**  
OB/GYN/Vascular Pkg. Color Doppler, PW Doppler, M-Mode, THI, DIMAQ, Image Data Management  
CD-Writer, Mitsubishi P93 B/W printer, DICOM, **CH5-2** Convex **EC9-4** Endo-cavity, **VF-13-5** Linear  
**w/ (1) year warranty -- for US / 220/110V -- ZONE C**
- 95 **Siemens CV-70 Cardio-Vascular (MFG: 2003) REV 4.0**  
Color Doppler, CW/PW Doppler, M-Mode, THI (harmonics), ECG/AUX, Dimaq IP, **STRESS ECHO**  
**P4-2** Cardiac Sector, **L10-5** Vascular/Small Parts Linear  
**110/220V -- ZONE A**  
"Special price" 3 units available at this price
- 96 **Siemens CV-70 Cardio-Vascular (MFG: 2003) REV 4.0**  
Color Doppler, CW/PW Doppler, M-Mode, THI (harmonics), ECG/AUX, Dimaq IP  
**P4-2** Cardiac Sector,  
**110/220V -- ZONE A**
- 97 **Siemens CV-70 Cardio-Vascular (MFG: 2003) REV 4.0**  
Color Doppler, CW/PW Doppler, M-Mode, THI (harmonics), ECG/AUX, Dimaq IP, **STRESS ECHO**,  
**DIMAQ IP, DICOM Print & Store**, Panasonic SVHS VCR, and Mitsubishi P91 B/W Printer  
**P4-2** Cardiac Sector,  
**110V/200V -- ZONE A**
- 98 **Siemens CV-70 Cardio-Vascular (MFG: 2003) REV 4.0**  
Color Doppler, CW/PW Doppler, M-Mode, THI (harmonics), **STRESS ECHO, DICOM, ECG/AUX**,  
**DIMAQ-IP, Panasonic SVHS VCR, Mitsubishi P91 B/W Printer**  
**P4-2** Cardiac Sector, **L10-5** Vascular/Small Parts Linear  
**110/220V -- ZONE A**
- 99 **Siemens CV-70 Cardio-Vascular (MFG: 2004) REV 4.0**  
Color Doppler, CW/PW Doppler, M-Mode, THI (harmonics), Dimaq IP, **STRESS ECHO**,  
**DICOM Print / Store / MWL**, ECG/AUX, Dicom MPPS, Axius-EF, **SieScape**, Image Arena M&R  
Panasonic SVHS VCR, Mitsubishi P91 B/W Printer  
**P4-2** Cardiac Sector, **L10-5** Vascular/Small Parts Linear  
**110/220V -- ZONE A**
- 100 **Siemens CV-70 Cardio-Vascular (MFG: 2005) REV 4.0**  
Color Doppler, CW/PW Doppler, M-Mode, THI (harmonics), Dimaq IP, **STRESS ECHO**,  
**DICOM Print / Store / MWL**, ECG/AUX, Dicom MPPS, Axius-EF, **SieScape**, Image Arena M&R  
Panasonic SVHS VCR, Mitsubishi P91 B/W Printer  
**P4-2** Cardiac Sector, **L10-5** Vascular/Small Parts Linear  
**ZONE A**

## **MINDRAY**

- 101 **MindRay DP-6600 (Portable) – MFG: 2006 – Demo Unit**  
B/W system, USB Port for Image Storage on to Memory Stick, Dual Probe Port,  
10" Monitor, **35C50EA - 3.5MHz Convex/50mmR – NEW PROBE**  
**ZONE A 1 year warranty**
- 102 **MindRay DP-9900 PLUS (Stand-Alone) – Demo Unit – MFG: 2007**  
B/W System, USB Port for Image Archive, CD Writer, 12" Monitor, Dual Probe Port, Swivel Keyboard  
Tissue Harmonic Imaging, Multi-Frequency, 256 Frame Cine, Built-in Image Memory  
**35C50HA (2.5/3.5/5.0MHz) Convex, 65EC10HA (5.0/6.5/7.5) MHz Endo-cavity**  
**ZONE A 1 year warranty**

**KPI Ultrasound, Inc.**  
**887 Marlborough Ave**  
**Riverside, CA 92507**

**Tel: (951) 367-0872**  
**Fax: (951) 367-0891**  
[www.kpiultrasound.com](http://www.kpiultrasound.com)

- 103 **MindRay DC-6 Color Doppler – ON/GYN/Vascular Calc. Pkg. --- NEW Unit**  
15" Monitor, Color, PW Doppler, Power Doppler, Directional Power Doppler  
Tissue Harmonics, Built-In Image Archive, CD-R/W and USB Ports  
**3C5** (2.5-5.0MHz) Convex, **7L4** (5.0-10.0MHz) Linear, **6CV1** (5.0/6.5/8.0MHz) Endo-cavity  
**SOLD** 1 year warranty

### VCR's (Reconditioned)

- **Panasonic AG-MD-830 SVHS VCR:**
- **Panasonic AG-MD-835 SVHS VCR:**
- **Sony SVO-9500MD SVHS VCR:**
- **Mitsubishi 3000MD SVHS VCR:**

Special: \*Mitsubishi MD-3000 SVHS VCR

### Color Printers (Reconditioned)

- **Sony UP-5200:**
- **Sony UP-2900:**
- **Sony UP-5600:**
- **Sony UP-51MD:**
- **Mitsubishi CP-700:**
- **Mitsubishi CP-800:**

### B/W Printers (Reconditioned)

- **Sony UP-910MD (Large Format):**
- **Sony UP-960MD (large Format):**
- **Sony UP-870MD**
- **Sony UP-890MD**
- **Sony UP895MD**
- **Sony UP-897MD (new)**
- **Sony UPD-895MD (new), Digital**
- **Sony UPD-897MD (new), Digital**
- **Mitsubishi P91**
- **Mitsubishi P93**

### MEDI-CAPTURE (Digital Archiving Unit)

- **USB-150 (still frame with LCD display)**
- **USB-200**  
(AVI, Still Frame, Movie with LCD display, Dicom)